

DATASHEET

The Baby Buddy manikins are lightweight, portable and suitable for an economical and realistic infant CPR training.

They feature a realistic head tilt and chest rise, palpable landmarks, a sanitary airway system and a lightweight and compact design.

Sanitary system

The disposable lung/mouth protection bags are changed with a lung insertion tool after each student preventing cross-contamination and the need to sanitize the manikins.

Easy AED training

The closed and anatomically correct torsi allow the placement of defibrillation pads for realistic, simulated AED training.

Note: Use only AED pads designed for foam body manikins.

The Baby Buddy features a removable chest for an easy installation of the protection lung bag.

Baby Buddy®



KEY BENEFITS

- Practice infant CPR techniques
- Lightweight, portable and economical
- Tough and durable
- Anatomically correct with palpable landmarks
- Airways open only with correct hyperextension
- Visible chest rise during ventilation
- Rebreathing airway
- Sanitary mouth/lung protection bags
- Easy AED training

ORDERING SPECIFICATIONS

Item no.

LF03721U

Baby Buddy 5-Pack

5 infant full-body manikins, 50 lung/mouth protection bags, 5 lung insertion tools, carry bag

LF03723U

Face Shield/Lung System for Sani-Child (set of 100)

Ambu

📍 Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com